


Nomen : _____

Datum : _____



_____ mm Diameter pupillae

Relat. Intens.

4					3					2					1					
e d c b a					e d c b a					e d c b a					e d c b a					
dB	0	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19
0																				
I																				
II																				
III																				
IV																				
V																				

Object

Correctio: _____ °

Refractio: _____ sph C _____ cyl _____ Visus: _____

OS. OD.

Für Seitenwechsel
Zeiger hier durchführen

Pour changer de côté,
faire passer l'index par ici

To change the side,
swing index along this line

Relat. Intens.	dB	Object
1	0,0315	15
2	0,100	10
3	0,315	5
4	1,00	0
a	0,40	I
b	0,50	II
c	0,63	III
d	0,80	IV
e	1,00	V

mm²

1

4

16

64